

Application Plate (AP-P)

For the efficient loading of IQ Ultra™ and IQ Ultimate™ Patch Test Units

- The device significantly facilitates preloading of multiple **IQ Ultra™** or **IQ Ultimate™** Patch Test Units. This is convenient in preparing test series in advance, such as the Baseline series, for a suitable number of patients.

- The improved design prevents the tape unit from sliding during application of the haptens.

- The device is equipped with a special function to detach and attach the syringe caps. This eliminates contamination of fingers during the dispensing of haptens.

- The device has a low weight of only 0.3 kg.

- The dimensions are (w x h): 43 x 27 cm.

