

Chemotechnique Patch Test Products - Amendments February 2023

New Series

Art. No.	Compound / Description	Series / Notes
PP-1000	Photopatch Series	Complete Series

Amended Haptens

Art. No.	Compound / Description	Series / Notes
A010	Ammonium hexachloroplatinate(IV) 0.1% aq ☒	Added to IMP-1000 Series - No. 25 (replaced A013) Added to MET-1000 Series - No. 19 (replaced A013) Deleted from METE-1000 Series (replaced with Z009)
A028	Acid Red 359 5% pet	Moved in TF-1000 Series - Moved to No. 5 (from No.32)
B003B	BENZISOTHIAZOLINONE (1,2-Benzisothiazolin-3-one, BIT) 0.1% pet ^{GB}	Added to S-1000 series - No. 30 (replaced MX30) Moved in ECB-1000 Series - moved to No. 30 (from No. 40)
B015B	2-BROMO-2-NITROPROPANE-1,3-DIOL (Bronopol) 0.5% pet ^{GB}	Moved in ECB-1000 Series - Moved to No. 33 (from No. 32)
C007A	QUATERNIUM-15 (Dowicil 200) 1% pet ^{GB}	Removed from S-1000 series (replaced with S011) Removed from ECB-1000 series (replaced with S011)
D044A	DIAZOLIDINYL UREA (Germall II) 2% pet ^{GB}	Moved in ECB-1000 Series - moved to No. 34 (from No. 33)
D065	DECYL GLUCOSIDE 5% pet ☒ ^{GB} **	Added to S-1000 Series - No. 32 Moved in ECB-1000 Series - to No. 32 (from No. 42)
H022	Hexamethylene diisocyanate (HDI) 0.1% pet ☒	Moved in I-1000 Series - moved to No. 2 (from No.6)
H031A	Hydroperoxides of Linalool 1% pet ☒ ^{GB} ****	Moved in ECB-1000 Series - moved to No. 37 (from No. 36)
H031B	Hydroperoxides of Linalool 0.5% pet ☒ ^{GB} ****	Moved in ECB-1000 Series - moved to No. 38 (from No. 37)
H032A	Hydroperoxides of Limonene 0.3% pet ☒ ^{GB} ****	Moved in ECB-1000 Series - moved to No. 39 (from No. 38)
H032B	Hydroperoxides of Limonene 0.2% pet ☒ ^{GB} ****	Moved in ECB-1000 Series - moved to No. 40 (from No. 39)
I001A	IMIDAZOLIDINYL UREA (Germall 115) 2% pet ^{GB}	Deleted from ECB-1000 Series (replaced with D044A)
L004	LAURYL POLYGLUCOSE (Lauryl glycoside) 3% pet ☒ ^{GB}	Deleted from ECB-1000 series
MX14	Wood tar mix (pine, beech, juniper, birch) 12% pet ☒	Moved in LU-1000 Series - moved to No. 15 (from No. 27)
MX26	Disperse Blue mix 106/124 1% pet ^{GB} *	Moved in TF-1000 Series - moved to No. 24 (from No. 33)
MX29A	Compositae mix II 5% pet ☒ *	Added to ECB-1000 Series - No. 36 (replaced H031A)
MX29B	Compositae mix II 2.5% pet ☒ ^{GB} *	Deleted from ECB-1000 Series (replaced with O004)
MX30	Textile dye mix 6.6% pet *	moved in ECB-1000 Series - moved to No. 31 (from No. 30) Moved in TF-1000 Series - moved to No. 27 (from No. 34) Moved in S-1000 series - moved to No. 31 (from No. 30) No longer available as individual hapten (only as part of complete series)
O004	2-n-Octyl-4-isothiazolin-3-one 0.1% pet ☒ ^{GB}	Moved in ECB-1000 Series - moved to No. 35 (from No.41)
P001	Palladium(II)chloride 2% pet	Added to MET-10000 Series - No. 23 Deleted from METE Series (replaced by R013)
P042	PANTHENOL 5% pet	Moved in H-1000 Series - moved to No. 12 (from No.38)
R013	Rhodium(III)chloride hydrate 2% pet	Deleted from MET-1000 Series (replaced by P001) Added to METE-1000 Series - No. 19 (replaced P001)
S004	SORBITAN OLEATE (Span 80) 5% pet	Added to ECB-1000 Series - No. 42 (replaced moved D065)
S005	SORBITAN SESQUIOLEATE 20% pet	Added to ECB-1000 Series - No. 41 (replaced O004)
S011	SODIUM METABISULFITE 1% pet ^{GB}	Added to S-1000 Series - No. 21 (replaces C007A) Moved in ECB-1000 Series - moved to No. 21 (from No. 31)
Z009	ZIRCONIUM DIOXIDE 0.1% pet	Moved in METE-1000 Series - moved to No. 18 (from No. 25)

Deleted Series

Art. No.	Compound / Description	Series / Notes
EP-1000	(Meth) Acrylate Series - Printing	Removed - some haptens from series available in other series
EPE-1000	(Meth) Acrylate Series - Printing	Removed - some haptens from series available in other series

Key to Notes: ☒ Hapten not recommended for pre-loading ^{GB} Present in BSCA Standard Series GB-1000 * Mix constituents available as individual haptens
** Emulsifier; Sorbitan Sesquioleate **** H031 & H032 have restricted shelf life & are maximum 3 hours out of the fridge ☼ Short dated & freezer storage

Deleted Haptens

Art. No.	Compound / Description	Series / Notes
A013	Ammonium tetrachloroplatinate(II) 0.25% aq ☒	Removed from PG-1000 series
B041	Benzydamine hydrochloride 2% pet	Removed from EP-1000 series
C004	Chlorhexidine diacetate 0.5% aq ☒	Removed from Various Haptens
C016	Coal tar 5% pet	Removed from Various Haptens
D005A	Dibucaine hydrochloride 5% pet	Removed from Various Haptens
D023B	Diphenylmethane-4,4'-diisocyanate (MDI) 0.5% pet ☒ **	Removed from IS-1000 series
E022	Eosin 5% pet	Removed from LU-1000 series
F006	Fenofibrate 10% pet	Removed from EPE-1000 series
I010B	Ibuprofen 5% pet	Removed from EPE-1000 series
O008	Olaquinox 1% pet	Removed from EPE-1000 series
P038	Polymeric diphenylmethane diisocyanate (PMDI) 2% pet ☒ **	Removed from I-1000 series

Amended Series

Art. No.	Compound / Description	Series / Notes
ECB-1000	European Comprehensive Baseline Series	No. 21 - S011 moved from No. 31 - replaces C007A (deleted from series) No.30 - B003B moved from No. 40 - replaces MX30 (moved to No. 31) No. 31 - MX30 moved from No. 30 - replaces S011 (moved to No. 21) No. 32 - D065 moved from No. 42 - replaces B015B (moved to No. 33) No. 33 - B015B moved from No. 32 - replaces D044A (moved to No. 34) No. 34 - D044A moved from No. 33 - replaces D044A (moved to No. 33) No. 35 - O004 moved from No. 41 - replaces O004 (deleted from series) No. 36 - MX29A added to series - replaces H031A (moved to No. 37) No. 37 - H031A moved from No. 36 - replaces H031B (moved to No. 38) No. 38 - H031B moved from No. 37 - replaces H032A (moved to No. 39) No. 39 - H032A moved from No. 38 - replaces H032B (moved to No. 40) No. 40 - H032B moved from 39 - replaces B003B (moved to No. 30) No. 41 - S005 added to series - replaces O004 (moved to No. 35) No. 42 - S004 added to series - replaces D065 (moved to No. 32) No. 43 - Empty (L004 deleted from series)
H-1000	Hairdressing Series	No. 12 - P042 moved from No. 38 - replaces empty space No. 38 - Empty (P042 moved to No. 12)
I-1000	Isocyanate series	No. 2 - H022 moved from No. 6 - replaces D023B (deleted from series) No. 6 - Empty (H022 moved to No. 2) No. 7 - Empty (P038 deleted from Series)
IMP-1000	Implant Series	No. 25 - A010 added to series (replaces A013 - deleted from series) No. 29 - P001 added to series (replaces R013 - deleted from Series)
LU-1000	Leg Ulcer Series	No. 15 - MX14 moved from No. 27 (replaces E022- deleted from series) No. 27 - Empty (MX14 moved to No. 15)
MET-1000	Metal Series	No. 19 - A010 added to series (replaces A013 - deleted from series) No. 23 - P001 added to series (replaces R013- deleted from Series)
METE-1000	Metal Extended Series	No. 18 - Z009 added to series (replaces A010 - deleted from series) No. 19 - R013 added to series (replaces P001 - deleted from series) No. 25 - Empty (Z009 moved to No. 18)
S-1000	European Baseline Series	No. 21 - S011 added to series (replaces C007A - deleted from series) No. 30 - B003B added to Series (replaces MX30 - moved to No. 31) No. 31 - MX30 moved from No. 30 No. 32 - D065 added to series
TF-1000	Textile Colours & Finish	No. 5 - A028 moved from No. 32 (replaces empty space) No. 24 - MX26 moved from No. 33 (replaces empty space) No. 27 - MX30 moved from No. 34(replaces empty space) No. 32 - Empty (A028 moved to No. 5) No. 33 - Empty (A028 moved to No. 24) No. 34 - Empty - A028 moved to No. 27
V-1000	Various Haptens	No. 3 - C004 deleted from series No. 5 - D005A deleted from series

Key to Notes: ☒ Hapten not recommended for pre-loading **GB** Present in BSCA Standard Series GB-1000 * Mix constituents available as individual haptens
 ** Emulsifier; Sorbitan Sesquioleate **** H031 & H032 have restricted shelf life & are maximum 3 hours out of the fridge ✱ Short dated & freezer storage